

Chemical Composition of Ionogem

Ionogem Powder ;

Alumino silicate fluoride barium glass	92.995%
Yellow Iron oxide	0.035%
Calcium tungstate	7%

Ionogem Liquid ;

Polyacrylic Acid	35%
Water	65%

Powder; component

	% w/w
Vitrified mixture of calcium/sodium fluorides Aluminium oxides, Phosphates and silicates	99.9
Iron oxide pigments	0.1

Liquid; Polyacrylic acid

Water	65.0
-------	------

Product code; 1180