

Chemical Composition of Gem Base LC

Gem Base Powder

Polyacrylic	15.10%
Barium Alumino silico phospho- fluoride glass	75.57%
Stabiliser	1.97%
Light Initiator	5.99%
Pigments	1.47%

Gem Base Liquids

Water	18.75%
Arylic Acid	16.52%
Hydroxyethyl Methacrylate	35.08%
Hexanedioldi Methacrylate	28.47%
Light Catalyst	1.18%

Product code ; 2140

Product code ; 2110