

## Chemical Composition of Gem Seal LC

### Formulation

	w/w
Butylated Hydroxy Toulene	0.02%
Bisphenol A Glyidyl Dimethacrylate	50.55%
Urethane Dimethacrylate	28.90%
Camphoroquinone	0.93%
Dimethyl Amino Ethyl Methacrylate	0.81%
Ethanol	9.00%
Tetrahydrofurfuryl-methacrylate	4.07%
Fumed silica	3.57%
Pigments ( EEC approved)	2.15%

Product code; 2200

Product code ; 2110