

Chemical Composition of Gem Line

Activator Paste

| | w/w |
|---|--------|
| Hexandiol 1-6 Dimethacrylate | 8.97% |
| Tetrahydrofurfuryl 2-Methacrylate | 2.85% |
| Dimethacrylate Amino Ethyl Methacrylate | 0.81% |
| Camphoroquinone | 0.94% |
| Zinc oxide | 0.53% |
| Pigments (EEC approved) | 0.09% |
| Fluridated glass | 48.65% |
| Bisphenol A Glycidyl Dimethacrylate | 33.48% |
| Fumed Silica | 3.68% |

Base paste

| | |
|-------------------------------------|--------|
| Polyacrylic Acid | 31.13% |
| Pigment | 3.00% |
| Butylated Hydroxy toluene | 0.06% |
| Benzoyl Peroxide | 0.60% |
| Bisphenol A Glycidyl Dimethacrylate | 63.21% |
| Fumed silica | 2.00% |

Product code ; 2120

Product code ; 2110