

Chemical Composition of Gem Core - Blue

Gem Core Powder

Alumino silicate fluoride barium glass	92.9955%
Blue iron oxide	0.0045%
Calcium tungstate	7%

Gem Core Liquid

Polyacrylic Acid	35%
Water	65%

Ferric Ferrocyanide C I (1971) No. 77510

Pigment Blue 27 complies with FDA Specifications

Above is similar in indications to ESPE keta silver and chemfil (detray, West Germany)

Product ; Ceramcore

Product code; 1130