

## Chemical Composition of Ortho Bond

<b>Paste</b>	w/w
Benzoyl Peroxide	4.90%
Butylated Hydroxy Toluene	0.14%
Hexane 1-6-Diol Dimethacrylate	18.90%
Bisphenol A Glycidyl Dimethacrylate	31.46%
Silica	44.60%

<b>Primer</b>	
Urethane Dimethacrylate	68.96%
Bisphenol A Glycidyl Dimethacrylate	19.08%
Tetra Hydro Furfuryl Methacrylate	11.10%
UV stabiliser	0.58%
Butylated Hydroxy Toluene	0.08%

Product code; 3100

Product code ; 2110