

Chemical Composition of Pulp Plus

Catalyst Paste

| | w/w |
|-----------------------|-------|
| Polymeric plasticiser | 42.0% |
| Zinc oxide | 20.0% |
| Silica | <1.0% |
| Calcium hydroxide | 37.0% |
| Titanium dioxide | <1.0% |

Base Paste

| | |
|--------------------|-------|
| Salicylate ester | 42.0% |
| Titanium dioxide | 4.0% |
| Calcium tungstate | 7.0% |
| Calcium phosphate | 31.0% |
| Caster oil | 7.0% |
| Barium sulphate | 7.0% |
| Iron oxide pigment | <1.0% |
| Silica | <1.0% |

Product code ; 2110