

## **PHYSICAL DATA**

### **ETCHING GEL/LIQUID**

#### **Description**

Dentine/Enamel Phosphoric Acid Etchant

#### **Active Ingredients**

38 - 38.22% w/w O-Phosphoric Acid

#### **Dose ratio**

Soaked cotton and pladget

#### **Etching time**

Permanent teeth Enamel 90 seconds

Deciduous teeth enamel 60 seconds

#### **Appearance**

Blue Thixotropic Gel/liquid

Liquid Viscosity 5000 - 6500 cp

Gel Viscosity 11000 - 12500 cp

Products code ; 5120