

PHYSICAL DATA

GEM BASE LC

FEATURES

Type 111 VLC radiopaque glass ionomer base/liner with superior properties due to presence of select acrylic and methacrylic acids.

Presented as a powder/liquid formation it gives a smooth non sticky mix, ideal for spreading into thin uniform linings or for base filling cavities prior to composite placement.

Ceramlux B/L is less moisture sensitive than the conventional glass ionomer and further it can be directly placed into cavity without acid etching. When placing into very deep cavities, the placement of calcium hydroxide liner is recommended.

Cerramlux B/L rapidly cures under light to provide a strong base which can be acid etched, prior to coating with composite bond.

Shade of ceramlux is dentine yellow, and is radiopaque.

Also available in capsules.

product code; 2140