

PHYSICAL DATA

PRODUCT SPECIFICATION

PRODUCT; **GEM CEM**

STANDARD TESTED TO; ISO 9917 ; 1991 - Luting

TEST CONDITIONS

Environmental control of 23'C = 1'C Room Temperature

Relative Humidity 50% = 5%

REQUIREMENTS

Test Parameter	Pass Criteria
Net setting time - minutes	2.5 min 8.0 max
Compressive strength Mpa	70 min
Acid erosion. Mm/h	0.05 max
Opacity. C0.70	–
Acid-soluble As. Content mg/kg	2.0 max
Acid-soluble Pb content mg/++-0kg	100.0 max

Product code; 2160