

PHYSICAL DATA

GEM BOND

Material : Dentine / Enamel bonding resin

Photo - Cure and chemical cure

Adhesion ; Over 5mpa. Onto acid etched human enamel (limit ; 4. 0 Mpa minimum) .

Chemical Cure ; 20 - 40 seconds (Limit; less than 40 seconds)

Shelf Life ; 18 months at normal temperature
36 months at storage temperature of 5 - 12 °C

The above product is tested to the following specification; ASP/SM/11/94

product code; 2100

