

PHYSICAL DATA

PRODUCT SPECIFICATION

PRODUCT; VARNISH

Description

Dental varnish for moisture proofing restoration and glazing

Active Ingredients

12.0 % w/w natural gums

4.8 % w/w cellulose

Dose ratio

Soaked cotton wool padlet

Drying time

0.5 ml at 23C drying within 60 seconds down to a water proof film

Appearance

Clear to translucent yellow, low viscosity volatile solvent/resin mix.

Quick drying to a water proof film.

Prod code ; 5140