

## PHYSICAL DATA

## IONO GEM

Material ; Glass Ionomer Filling Material

Complies with ; BS 6039 ; 1981 Type II  
ISO 7489;1986 Type II

### Physical Property Data Test

Consistency	30 mm + 2mm
Film Thickness	-
Setting time	4 mins - 5 mins 50 sec
Working time	2 mins 15 sec - 2 mins 30 sec
Compressive Strength	143.24 Mpa
Diametral tensile strength	12.07 mpa
Water leachable cement	Pass
Acid soluble arsenic content	Pass

Products code; 1180