

## **PHYSICAL DATA**

### **PRODUCT SPECIFICATION**

**PRODUCT;** **CAVFIL PINK WHITE**

**STANDARD TESTED TO;** PSP In house standard  
Ref; 2602, 2614, 0692/ASP

**TEST CONDITIONS** Room temperature 23°C = 1°C  
Relative Humidity 50% = 5%

### **TEST PARAMETER**

#### **Physical Appearance**

The putty should have firm pliable consistency with no lumping of Powder or streaking of colour.

#### **The colour should be;**

For pink putty - between salmon pink and gum pink in shade  
For white putty- between pure white to light grey/yellow in shade

#### **Odour**

The putty should have a distinct mint odour

#### **Hardening time**

The hardening time is to be 2 minutes 45 seconds = 45s

#### **Cracking**

Only the smallest hair line crack is acceptable in not more than 1 of 3 test samples of the same batch .

cc QA 06/1, MDR Cavfil